



## F4/80

Biolegend (Cat #123102); or Novus (Cat # NBP2-12506) Immunohistochemistry Protocol For Formalin Fixed Frozen Embedded Tissue Cut sections at 7 microns and stored in -80°C

Day 1

- 1. Thaw the slides at room temperature in a humid chamber for 10 minutes.
- 2. Rinse the slides with 3 changes of PBS for 10 minutes each.
- 3. Outline the tissues with a **PAP pen**.
- 4. Rinse the slides with PBST (0.2% TritonX diluted in PBS) once for 10 minutes.
- 5. Incubate with **5% normal donkey serum** diluted in PBST (0.2% TritonX diluted in PBS) at room temperature in a humid chamber for **30 minutes**.
- 6. Rinse the slides with 3 changes of PBS for 10 minutes each.
- 7. Prepare primary antibody (1. *RatF4/80, Biolegend 123102*; or 2. *RabbitF4/80, Novus NBP2-12506*) in **1% normal donkey serum** diluted in PBS.
- 8. Incubate the slides with a **1:50** dilution of F4/80 primary antibody in a humid chamber at **4**°**C** overnight.

Day 2

- 9. Rinse the slides with 3 changes of PBS for 10 minutes each.
- 10. Prepare secondary antibody (1. *GoatαRat-Fluro488*, *Abcam ab150165*; or 2. *GoatαRabbit-Fluro488*, *Abcam ab150077*) in **1% normal donkey serum** diluted in PBS.
- 11. Incubate the slides with a **1:500** dilution of secondary antibody in a humid chamber at **room** temperature for 1 hour.
- 12. Rinse the slides with 3 changes of PBS for 10 minutes each.
- 13. Mount with **mounting medium containing Dapi** (Vectashield H-1200).

