

p13E-11 hybridizations**Dextrane Hybridization buffer (1 Liter):**

500 mL H₂O

200 mL 5M NaCl

25 mL 2M Tris-HCl (pH 7.5)

50 mL 20% SDS

100 g Dextran sodium sulfate (→add slowly, while mixing the solution)

Add until 1 Liter with demineralized water

Heat in water bath until 65°C, then stir until Dextran sulfate is dissolved

For hybridization:

Heat until 65°C

Add 100 ug/ml Salmon sperm DNA to hybridization buffer (preheated at 95°C) and start prehybridization (standard prehybridization at least 45 minutes)

We do not replace the 'pre'-hybridization buffer by fresh hybridization buffer.

D4Z4 hybridizations**2xNaPi hybridization buffer (1 liter),****use in combination with an equal volume of deionized Formamide)**

Prepare 2x phosphate buffer as follows:

approx. 500 mL demineralized water in 2 Liter beaker

Dissolve 44,45 gram Na₂HPO₄ * 2H₂O (MW=177,99, Merck, Fluka)

Dissolving takes 1 hour, heat a little bit, not more than 40°C

Set pH with 85% phosphoric acid (2-3 mL) to exactly pH = 7.2

Do not sterilize, prepare prior to preparation hybridization buffer.

Add to phosphate buffer:

100 mL NaCl

4 mL 0.5 M EDTA (pH 8.0)

140 gram SDS, or 700 mL 20% SDS solution (For solid SDS use the hood (very toxic))

Add until 1 Liter with demineralized water

Heat in water bath until 65°C, then stir until SDS is dissolved.

For hybridization:

Heat until 65°C

Add 1 volume of deionized Formamide

Add 100 ug/ml Salmon sperm DNA to hybridization buffer (preheated at 95°C) and start prehybridization (standard prehybridization at least 45 minutes)

We do not replace the 'pre'-hybridization buffer by fresh hybridization buffer.

For 4qA+4qB and all other hybridizations:

NaPi/PEG hybridization buffer (2 liter):

Prepare phosphate buffer (0.5 M) as follows:

approx. 500 mL demineralized water in 2 Liter beaker

Dissolve 44,45 gram $\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$ (MW=177,99, Merck, Fluka)

Dissolving takes 1 hour, heat a little bit, not more than 40°C

Set pH with 85% phosphoric acid (2-3 mL) to exactly pH = 7.2

Do not sterilize, prepare prior to preparation hybridization buffer.

Add to phosphate buffer:

100 mL NaCl (5M)

4 mL 0.5 M EDTA (pH 8.0)

140 gram SDS, or 700 mL 20% SDS solution (for solid SDS use the hood (very toxic))

200 gram PEG-6000

Add until almost 2 Liter with demineralized water

End concentration:

0.125 M phosphate buffer, 0,25M NaCl, 1mM EDTA, 7% SDS, 10% PEG6000

Heat in water bath until 65°C, then stir until PEG and SDS are dissolved.

Finally add demineralized water until 2 Liter.

For hybridization:

Heat until 65°C

Add 100 ug/ml Salmon sperm DNA to hybridization buffer (preheated at 95°C) and start prehybridization (standard prehybridization at least 45 minutes)

We do not replace the 'pre'-hybridization buffer by fresh hybridization buffer.