

## **Guidelines for Use of the BIAcore T100**

1. Please sign-up for use of the instrument, both for purposes scheduling and fees. The requested limit is one week block per laboratory reserved on the calendar at a time. As soon as the scheduled use is complete, that laboratory is welcome to reserve further use.

[http://www.brownbearsw.com/freecal/UofR\\_Biacore](http://www.brownbearsw.com/freecal/UofR_Biacore)

2. First time users of the U of R T100 instrument should be trained by Dr. Jenkins (Jermaine\_jenkins@urmc.rochester.edu) until he provides a 'certificate for independent use'. This training may be short (1-2 hours) or long (>1 day) depending on the experience of the user. **Dr. Jenkins will provide the password to the sign-up website at that time.**
3. Please run a 'desorb' cleaning cycle at the completion of instrument use, and return the instrument to standby mode with a large volume of water. Dr. Jenkins will run monthly sanitization procedures and replenish the supply of water for standby mode.

For more information, please see:

[http://dbb.urmc.rochester.edu/bcbp/facilities/bsb\\_facility.htm](http://dbb.urmc.rochester.edu/bcbp/facilities/bsb_facility.htm)